

For Research Use Only

NIACR1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83044-1-PBS

Basic Information

Catalog Number: 83044-1-PBS	GenBank Accession Number: BC027965	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 338442	CloneNo.: 230163B7
Source: Rabbit	UNIPROT ID: Q8TDS4	
Isotype: IgG	Full Name: niacin receptor 1	
Immunogen Catalog Number: AG12425	Calculated MW: 363 aa, 42 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00117-1, NIACR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83044-1-PBS. Detection antibody: 83044-4-PBS. Standard: Ag12425. Range: 0.625-20 ng/mL