

For Research Use Only

# CLDN5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83046-3-PBS

## Basic Information

Catalog Number: 83046-3-PBS	GenBank Accession Number: BC032363	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 7122	CloneNo.: 230250E5
Source: Rabbit	UNIPROT ID: O00501	
Isotype: IgG	Full Name: claudin 5	
Immunogen Catalog Number: AG29044	Calculated MW: 218 aa, 23 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00128-2, CLDN5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83046-1-PBS. Detection antibody: 83046-3-PBS. Standard: Ag29044. Range: 0.313-20 ng/mL