

For Research Use Only

KIAA0776 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83047-1-PBS

Basic Information

Catalog Number:

83047-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23362

GenBank Accession Number:

BC036379

GeneID (NCBI):

23376

UNIPROT ID:

O94874

Full Name:

KIAA0776

Purification Method:

Protein A purification

CloneNo.:

230363B4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00135-2, KIAA0776 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83047-1-PBS. Detection antibody: 83047-4-PBS. Standard: Ag23362. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00135-3, KIAA0776 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83047-1-PBS. Detection antibody: 83047-2-PBS. Standard: Ag23362. Range: 0.625-20 ng/mL