

For Research Use Only

KIAA0776 Recombinant antibody, PBS Only (Capture)



Catalog Number:83047-3-PBS

Basic Information

| | | |
|--------------------------------------|---------------------------------------|--|
| Catalog Number: 83047-3-PBS | GenBank Accession Number: BC036379 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 23376 | CloneNo.: 230363E9 |
| Source: Rabbit | UNIPROT ID: O94874 | |
| Isotype: IgG | Full Name: KIAA0776 | |
| Immunogen Catalog Number: AG23362 | | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00135-1, KIAA0776 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83047-3-PBS. Detection antibody: 83047-4-PBS. Standard: Ag23362. Range: 0.313-20 ng/mL.