

For Research Use Only

# KIAA0776 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83047-4-PBS

## Basic Information

Catalog Number: 83047-4-PBS	GenBank Accession Number: BC036379	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 23376	CloneNo.: 230363F10
Source: Rabbit	UNIPROT ID: O94874	
Isotype: IgG	Full Name: KIAA0776	
Immunogen Catalog Number: AG23362		

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00135-1, KIAA0776 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83047-3-PBS. Detection antibody: 83047-4-PBS. Standard: Ag23362. Range: 0.313-20 ng/mL.



Cytometric bead array standard curve of MP00135-2, KIAA0776 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83047-1-PBS. Detection antibody: 83047-4-PBS. Standard: Ag23362. Range: 0.313-20 ng/mL.