

For Research Use Only

Mouse EPO Recombinant antibody, PBS Only (Detector)



Catalog Number:83050-1-PBS

Basic Information

Catalog Number:

83050-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_007942

GeneID (NCBI):

13856

UNIPROT ID:

P07321

Full Name:

erythropoietin

Calculated MW:

21kd

Purification Method:

Protein A purification

CloneNo.:

230304A12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00026-1, Mouse EPO Recombinant Matched Antibody Pair, PBS Only. 83050-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0265. 83050-1-PBS was HRP conjugated as the detection antibody. Range: 7.8-500 pg/mL



Sandwich ELISA standard curve of MP00026-2, Mouse EPO Recombinant Matched Antibody Pair, PBS Only. 83050-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0265. 83050-1-PBS was HRP conjugated as the detection antibody. Range: 7.8-500 pg/mL