

For Research Use Only

# Integrin beta 3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83053-2-PBS

## Basic Information

<b>Catalog Number:</b> 83053-2-PBS	<b>GenBank Accession Number:</b> BC127667	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 3690	<b>CloneNo.:</b> 230292B6
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P05106	
<b>Isotype:</b> IgG	<b>Full Name:</b> integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61)	
	<b>Calculated MW:</b> 87 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00047-1, Integrin beta-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83053-2-PBS. Detection antibody: 83053-3-PBS. Standard: SY00181. Range: 0.312-5 ng/mL.



Cytometric bead array standard curve of MP00047-2, Integrin beta-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83053-2-PBS. Detection antibody: 83053-1-PBS. Standard: SY00181. Range: 0.312-5 ng/mL.