

For Research Use Only

Integrin beta 3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83053-3-PBS

Basic Information

Catalog Number: 83053-3-PBS	GenBank Accession Number: BC127667	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3690	CloneNo.: 230292F12
Source: Rabbit	UNIPROT ID: P05106	
Isotype: IgG	Full Name: integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61)	
	Calculated MW: 87 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00047-1, Integrin beta-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83053-2-PBS. Detection antibody: 83053-3-PBS. Standard: SY00181. Range: 0.312-5 ng/mL.