

For Research Use Only

OATP1B1/SLCO1B1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83054-4-PBS



Basic Information

Catalog Number:

83054-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_006446

GeneID (NCBI):

10599

UNIPROT ID:

Q9Y6L6

Full Name:

solute carrier organic anion
transporter family, member 1B1

Calculated MW:

76 kDa

Purification Method:

Protein A purification

CloneNo.:

230090H6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00051-2, SLCO1B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83054-4-PBS. Detection antibody: 83054-3-PBS. Standard: SY00964. Range: 0.312-10 ng/mL.