

For Research Use Only

CD14 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83055-3-PBS

Basic Information

Catalog Number: 83055-3-PBS	GenBank Accession Number: BC010507	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 929	CloneNo.: 230332E5
Source: Rabbit	ENSEMBL Gene ID: ENSG00000170458	
Isotype: IgG	UNIPROT ID: P08571	
	Full Name: CD14 molecule	
	Calculated MW: 375 aa, 40 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00058-3, CD14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83055-3-PBS. Detection antibody: 83055-5-PBS. Standard: Eg0498. Range: 0.625-20 ng/mL.