

For Research Use Only

EPOR Recombinant antibody, PBS Only (Detector)



Catalog Number: 83056-3-PBS

Basic Information

Catalog Number: 83056-3-PBS	GenBank Accession Number: NM_000121.4	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2057	CloneNo.: 230117F9
Source: Rabbit	UNIPROT ID: P19235	
Isotype: IgG	Full Name: erythropoietin receptor	
Immunogen Catalog Number: AG33956	Calculated MW: 55 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00065-1, EPOR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83056-4-PBS. Detection antibody: 83056-3-PBS. Standard: Ag33956. Range: 0.625-40 ng/mL