

For Research Use Only

GABRA2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83057-1-PBS

Basic Information

Catalog Number: 83057-1-PBS	GenBank Accession Number: NM_000807	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2555	CloneNo.: 230158C12
Source: Rabbit	UNIPROT ID: P47869	
Isotype: IgG	Full Name: gamma-aminobutyric acid (GABA) A receptor, alpha 2	
Immunogen Catalog Number: AG29301	Calculated MW: 51 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00067-1, GABRA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83057-3-PBS. Detection antibody: 83057-1-PBS. Standard: Ag29301. Range: 0.312-20 ng/mL.