

For Research Use Only

Cytokeratin 13 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83058-3-PBS

Basic Information

| | | |
|--------------------------------------------|----------------------------------------------|-------------------------------------------------------|
| Catalog Number: 83058-3-PBS | GenBank Accession Number: BC002661 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 3860 | CloneNo.: 230366F3 |
| Source: Rabbit | UNIPROT ID: P13646 | |
| Isotype: IgG | Full Name: keratin 13 | |
| Immunogen Catalog Number: AG0217 | Calculated MW: 50 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00074-1, KRT13 Recombinant Matched Antibody Pair - BSA and Azide free. Capture antibody: 83058-1-PBS. Detection antibody: 83058-3-PBS. Standard: Ag0217. Range: 0.625-20 ng/mL.