

For Research Use Only

SOX2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83059-1-PBS

Basic Information

Catalog Number:

83059-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0612

GenBank Accession Number:

BC013923

GeneID (NCBI):

6657

UNIPROT ID:

P48431

Full Name:

SRY (sex determining region Y)-box 2

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

230225A7

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00077-1 measured by cytometric bead array. Capture antibody: 83059-1-PBS. Detection antibody: 83059-3-PBS. Range: 0.625-20 ng/mL