

For Research Use Only

SOX2 Recombinant antibody, PBS Only (Detector)



Catalog Number:83059-3-PBS

Basic Information

Catalog Number: 83059-3-PBS	GenBank Accession Number: BC013923	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 6657	CloneNo.: 230225E7
Source: Rabbit	UNIPROT ID: P48431	
Isotype: IgG	Full Name: SRV (sex determining region Y)-box 2	
Immunogen Catalog Number: AG0612	Calculated MW: 34 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00077-1, SOX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83059-1-PBS. Detection antibody: 83059-3-PBS. Standard: Ag0612. Range: 0.625-20 ng/mL