

For Research Use Only

ABCA13 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83060-2-PBS

Basic Information

Catalog Number: 83060-2-PBS	GenBank Accession Number: NM_152701	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 154664	CloneNo.: 230288G4
Source: Rabbit	UNIPROT ID: Q86UQ4	
Isotype: IgG	Full Name: ATP-binding cassette, sub-family A (ABC1), member 13	
	Calculated MW: 576 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00084-1, ABCA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83060-1-PBS. Detection antibody: 83060-2-PBS. Standard: SY00509. Range: 0.312-20 ng/mL



Cytometric bead array standard curve of MP00084-2, ABCA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83060-3-PBS. Detection antibody: 83060-2-PBS. Standard: SY00509. Range: 0.312-20 ng/mL