

For Research Use Only

MOG Recombinant antibody, PBS Only (Capture)



Catalog Number: 83063-1-PBS

Basic Information

Catalog Number: 83063-1-PBS	GenBank Accession Number: BC035938	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4340	CloneNo.: 230291C10
Source: Rabbit	UNIPROT ID: Q16653	
Isotype: IgG	Full Name: myelin oligodendrocyte glycoprotein	
Immunogen Catalog Number: AG3273	Calculated MW: 295 aa, 34 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00110-1, MOG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83063-1-PBS. Detection antibody: 83063-4-PBS. Standard: Ag3273. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00110-3, MOG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83063-1-PBS. Detection antibody: 83063-3-PBS. Standard: Ag3273. Range: 0.625-20 ng/mL