

For Research Use Only

Mouse KIM-1/HAVCR1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83064-1-PBS

Basic Information

Catalog Number: 83064-1-PBS	GenBank Accession Number: NM_134248	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 171283	CloneNo.: 230336A5
Source: Rabbit	UNIPROT ID: Q5QNS5	
Isotype: IgG	Full Name: hepatitis A virus cellular receptor 1	
	Calculated MW: 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00116-2, Mouse KIM-1/HAVCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83064-1-PBS. Detection antibody: 83064-2-PBS. Standard: Eg0571. Range: 0.625-20 ng/mL.