

For Research Use Only

Mouse KIM-1/HAVCR1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83064-2-PBS

Basic Information

Catalog Number:

83064-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_134248

GeneID (NCBI):

171283

UNIPROT ID:

Q5QNS5

Full Name:

hepatitis A virus cellular receptor 1

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

230336B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00116-2, Mouse KIM-1/HAVCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83064-1-PBS. Detection antibody: 83064-2-PBS. Standard: Eg0571. Range: 0.625-20 ng/mL.