

For Research Use Only

# ICAM-1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83069-1-PBS

## Basic Information

Catalog Number:

83069-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC015969

GeneID (NCBI):

3383

ENSEMBL Gene ID:

ENSG00000090339

UNIPROT ID:

P05362

Full Name:

intercellular adhesion molecule 1

Calculated MW:

90 kDa

Purification Method:

Protein A purification

CloneNo.:

230342F1

## Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data



Standard curve of MP00028-5 measured by sandwich ELISA. Capture antibody: 83069-1-PBS. Detection antibody: 83069-2-PBS. Range: 39-2500 pg/mL



Standard curve of MP00028-1 measured by cytometric bead array. Capture antibody: 83069-1-PBS. Detection antibody: 83069-3-PBS. Range: 0.625-20 ng/mL



Standard curve of MP00028-2 measured by cytometric bead array. Capture antibody: 83069-1-PBS. Detection antibody: 83069-4-PBS. Range: 0.625-20 ng/mL