

For Research Use Only

ICAM-1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83069-2-PBS

Basic Information

Catalog Number:

83069-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC015969

GeneID (NCBI):

3383

ENSEMBL Gene ID:

ENSG00000090339

UNIPROT ID:

P05362

Full Name:

intercellular adhesion molecule 1

Calculated MW:

90 kDa

Purification Method:

Protein A purification

CloneNo.:

230342F2

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00028-4 measured by sandwich ELISA. Capture antibody: 83069-5-PBS. Detection antibody: 83069-2-PBS. Range: 78.1-5000 pg/mL



Standard curve of MP00028-5 measured by sandwich ELISA. Capture antibody: 83069-1-PBS. Detection antibody: 83069-2-PBS. Range: 39-2500 pg/mL