

For Research Use Only

# ICAM-1/CD54 Recombinant antibody, PBS Only (Detector)

Catalog Number:83069-3-PBS



## Basic Information

<b>Catalog Number:</b> 83069-3-PBS	<b>GenBank Accession Number:</b> BC015969	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 3383	<b>CloneNo.:</b> 230342F9
<b>Source:</b> Rabbit	<b>ENSEMBL Gene ID:</b> ENSG00000090339	
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P05362	
	<b>Full Name:</b> intercellular adhesion molecule 1	
	<b>Calculated MW:</b> 90 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00028-1, ICAM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83069-1-PBS. Detection antibody: 83069-3-PBS. Standard: Eg0570. Range: 0.625-20 ng/mL