

For Research Use Only

ICAM-1/CD54 Recombinant antibody, PBS Only (Detector)

Catalog Number:83069-4-PBS



Basic Information

| | | |
|---------------------------------------|--|---|
| Catalog Number: 83069-4-PBS | GenBank Accession Number: BC015969 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 3383 | CloneNo.: 230342H4 |
| Source: Rabbit | ENSEMBL Gene ID: ENSG00000090339 | |
| Isotype: IgG | UNIPROT ID: P05362 | |
| | Full Name: intercellular adhesion molecule 1 | |
| | Calculated MW: 90 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00028-2, ICAM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83069-1-PBS. Detection antibody: 83069-4-PBS. Standard: Eg0570. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00028-3, ICAM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83069-5-PBS. Detection antibody: 83069-4-PBS. Standard: Eg0570. Range: 0.625-20 ng/mL