

For Research Use Only

Mouse Adiponectin Recombinant antibody, PBS Only (Capture)



Catalog Number:83070-2-PBS

Basic Information

Catalog Number:

83070-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009605

GeneID (NCBI):

11450

UNIPROT ID:

Q60994

Full Name:

adiponectin, C1Q and collagen domain containing

Calculated MW:

27kd

Purification Method:

Protein A purification

CloneNo.:

230258C4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00039-1, Mouse Adiponectin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83070-2-PBS. Detection antibody: 83070-4-PBS. Standard: Eg0262. Range: 0.039-2.5 ng/mL