

For Research Use Only

Mouse Adiponectin Recombinant antibody, PBS Only (Detector)



Catalog Number:83070-4-PBS

Basic Information

| | | |
|---------------------------------------|--|---|
| Catalog Number: 83070-4-PBS | GenBank Accession Number: NM_009605 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 11450 | CloneNo.: 230258D8 |
| Source: Rabbit | UNIPROT ID: Q60994 | |
| Isotype: IgG | Full Name: adiponectin, C1Q and collagen domain containing | |
| | Calculated MW: 27kd | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00039-1, Mouse Adiponectin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83070-2-PBS. Detection antibody: 83070-4-PBS. Standard: Eg0262. Range: 0.039-2.5 ng/mL