

For Research Use Only

Mouse Adiponectin Recombinant antibody, PBS Only (Capture)



Catalog Number:83070-5-PBS

Basic Information

Catalog Number: 83070-5-PBS	GenBank Accession Number: NM_009605	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 11450	CloneNo.: 230258G11
Source: Rabbit	UNIPROT ID: Q60994	
Isotype: IgG	Full Name: adiponectin, C1Q and collagen domain containing	
	Calculated MW: 27kd	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00039-2, Mouse Adiponectin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83070-5-PBS. Detection antibody: 83070-1-PBS. Standard: Eg0262. Range: 0.039-2.5 ng/mL