

For Research Use Only

STT3A Recombinant antibody, PBS Only (Capture)



Catalog Number:83077-3-PBS

Basic Information

Catalog Number:

83077-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31994

GenBank Accession Number:

BC020965

GeneID (NCBI):

3703

UNIPROT ID:

P46977

Full Name:

STT3, subunit of the oligosaccharyltransferase complex, homolog A (*S. cerevisiae*)

Calculated MW:

705 aa, 81 kDa

Purification Method:

Protein A purification

CloneNo.:

230426G7

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00080-1 measured by cytometric bead array. Capture antibody: 83077-3-PBS. Detection antibody: 83077-2-PBS. Range: 0.625-20 ng/mL



Standard curve of MP00080-2 measured by cytometric bead array. Capture antibody: 83077-3-PBS. Detection antibody: 83077-1-PBS. Range: 0.625-20 ng/mL