

For Research Use Only

p300 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83078-2-PBS

Basic Information

Catalog Number: 83078-2-PBS	GenBank Accession Number: NM_001429	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2033	CloneNo.: 230369G4
Source: Rabbit	UNIPROT ID: Q09472	Affinity: $K_D = 2.37 \times 10^{-8}$
Isotype: IgG	Full Name: E1A binding protein p300	$K_{off} = 4.58 \times 10^{-3}$ $K_{on} = 1.93 \times 10^5$
	Calculated MW: 264 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C .
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00310-1, EP300 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83078-2-PBS. Detection antibody: 83078-1-PBS. Standard: SY00911. Range: 0.313-10 ng/mL