

For Research Use Only

STK39 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83081-2-PBS

Basic Information

Catalog Number:

83081-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34503

GenBank Accession Number:

BC166614

GeneID (NCBI):

27347

UNIPROT ID:

Q9UEW8

Full Name:

serine threonine kinase 39
(STE20/SPS1 homolog, yeast)

Calculated MW:

59 kDa

Purification Method:

Protein A purification

CloneNo.:

230382G9

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00106-2 measured by cytometric bead array. Capture antibody: 83081-2-PBS. Detection antibody: 83081-4-PBS. Range: 0.625-20 ng/mL