

For Research Use Only

STK39 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83081-4-PBS

Basic Information

Catalog Number: 83081-4-PBS	GenBank Accession Number: BC166614	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 27347	CloneNo.: 230382H4
Source: Rabbit	UNIPROT ID: Q9UEW8	
Isotype: IgG	Full Name: serine threonine kinase 39 (STE20/SPS1 homolog, yeast)	
Immunogen Catalog Number: AG34503	Calculated MW: 59 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00106-2, STK39 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83081-2-PBS. Detection antibody: 83081-4-PBS. Standard: Ag34503. Range: 0.625-20 ng/mL