

For Research Use Only

DDX39 Recombinant antibody, PBS Only (Capture)

Catalog Number:83083-1-PBS



Basic Information

Catalog Number:

83083-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2311

GenBank Accession Number:

BC032128

GeneID (NCBI):

10212

UNIPROT ID:

O00148

Full Name:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 39

Calculated MW:

427 aa, 49 kDa

Purification Method:

Protein A purification

CloneNo.:

230375B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00111-2, DDX39 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83083-1-PBS. Detection antibody: 83083-2-PBS. Standard: Ag2311. Range: 1.25-20 ng/mL