

For Research Use Only

DDX39 Recombinant antibody, PBS Only (Detector)

Catalog Number:83083-2-PBS



Basic Information

Catalog Number: 83083-2-PBS	GenBank Accession Number: BC032128	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 10212	CloneNo.: 230375E8
Source: Rabbit	UNIPROT ID: O00148	
Isotype: IgG	Full Name: DEAD (Asp-Glu-Ala-Asp) box polypeptide 39	
Immunogen Catalog Number: AG2311	Calculated MW: 427 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00111-1, DDX39 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83083-3-PBS. Detection antibody: 83083-2-PBS. Standard: Ag2311. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00111-2, DDX39 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83083-1-PBS. Detection antibody: 83083-2-PBS. Standard: Ag2311. Range: 1.25-20 ng/mL