### For Research Use Only

# DDX39 Recombinant antibody, PBS Only (Capture)



Catalog Number:83083-3-PBS

#### **Basic Information**

Catalog Number: 83083-3-PBS Size: 1mg/ml Source: Rabbit Isotype: IgG Immunogen Catalog Number: AG2311 GenBank Accession Number: BC032128 GeneID (NCBI): 10212 UNIPROT ID: 000148 Full Name: DEAD (Asp-Glu-Ala-Asp) box polypeptide 39 Calculated MW: 427 aa, 49 kDa Purification Method: Protein A purification CloneNo.: 230375G6

**Applications** 

Tested Applications: Indirect ELISA, Cytometric bead array Species Specificity: Human

#### **Background Information**


Storage

Storage: Store at -80°C. The product is shipped with ice packs. Upon receipt, store it immediately at -80°C Storage Buffer: 100% PBS pH 7.3

For technical support and original validation data for this product please contact:T: 4006900926E: Proteintech-CN@ptglab.comW: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00111-1, DDX39 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83083-3-PBS. Detection antibody: 83083-2-PBS. Standard: Ag2311. Range: 0.625-20 ng/mL