

For Research Use Only

# MOGAT2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83086-2-PBS

## Basic Information

Catalog Number:

83086-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15806

GenBank Accession Number:

BC103878

GeneID (NCBI):

80168

UNIPROT ID:

Q3SYC2

Full Name:

monoacylglycerol O-acyltransferase

2

Calculated MW:

334 aa, 38 kDa

Purification Method:

Protein A purification

CloneNo.:

230365F1

## Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00124-1 measured by cytometric bead array. Capture antibody: 83086-1-PBS. Detection antibody: 83086-2-PBS. Range: 0.625-20 ng/mL