

For Research Use Only

# GDNF Recombinant antibody, PBS Only



Catalog Number: 83087-2-PBS

## Basic Information

Catalog Number: 83087-2-PBS	GenBank Accession Number: BC069369	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2668	CloneNo.: 230358G6
Source: Rabbit	UNIPROT ID: P39905	
Isotype: IgG	Full Name: glial cell derived neurotrophic factor	
Immunogen Catalog Number: AG24068	Calculated MW: 211 aa, 24 kDa	

## Applications

Tested Applications:  
WB, ELISA  
Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 83087-2-RR (GDNF antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83087-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83087-2-RR against Human GDNF were performed. The affinity constant is below 1 pM.