

For Research Use Only

PCP4L1 Recombinant antibody, PBS Only



Catalog Number: 83093-2-PBS

Basic Information

Catalog Number: 83093-2-PBS	GenBank Accession Number: BC028905	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 654790	CloneNo.: 230203G12
Source: Rabbit	UNIPROT ID: A6NKN8	
Isotype: IgG	Full Name: Purkinje cell protein 4 like 1	
Immunogen Catalog Number: AG31230	Calculated MW: 68 aa, 7 kDa	

Applications

Tested Applications:
IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

PCP4L1 (Purkinje cell protein 4-like 1) is a small neuronal IQ motif protein closely related to the calmodulin-binding protein Pcp4/PEP-19 (PMID: 22458599).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed Caco-2 cells using PCP4L1 antibody (83093-2-RR, Clone: 230203G12) at dilution of 1:250 and CoraLite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 83093-2-PBS in a different storage buffer formulation.



1x10⁶ Caco-2 cells were intracellularly stained with 0.25 ug Anti-Human PCP4L1 (83093-2-RR, Clone:230203G12) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Rabbit IgG control Rabbit PolyAb (30000-0-AP) (blue). Cells were fixed and permeabilized with True-Nuclear Transcription Factor Buffer Set. This data was developed using the same antibody clone with 83093-2-PBS in a different storage buffer formulation.