

For Research Use Only

TGOLN2,TGN46 Recombinant antibody, PBS Only (Capture)



Catalog Number:83099-5-PBS

Basic Information

Catalog Number: 83099-5-PBS	GenBank Accession Number: BC028219	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 10618	CloneNo.: 230376F9
Source: Rabbit	UNIPROT ID: O43493	Affinity: $K_D=1.11 \times 10^{-9}$
Isotype: IgG	Full Name: trans-golgi network protein 2	$K_{off}=2.89 \times 10^{-4}$
Immunogen Catalog Number: AG4470	Calculated MW: 447 aa, 47 kDa	$K_{on}=2.59 \times 10^5$

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00152-2, TGOLN2,TGN46 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83099-5-PBS. Detection antibody: 83099-3-PBS. Standard: Ag4470. Range: 0.625-20 ng/mL.