

For Research Use Only

PCYT1A Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83100-4-PBS

Basic Information

Catalog Number: 83100-4-PBS	GenBank Accession Number: NM_005017	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5130	CloneNo.: 230393D8
Source: Rabbit	UNIPROT ID: P49585	
Isotype: IgG	Full Name: phosphate cytidyltransferase 1, choline, alpha	
Immunogen Catalog Number: AG34268	Calculated MW: 42 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00303-1, PCYT1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83100-3-PBS. Detection antibody: 83100-4-PBS. Standard: Ag34268. Range: 1.25-20 ng/mL.



Cytometric bead array standard curve of MP00303-2, PCYT1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83100-4-PBS. Detection antibody: 83100-2-PBS. Standard: Ag34268. Range: 1.25-20 ng/mL.