

For Research Use Only

# TNFR2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83101-4-PBS



## Basic Information

Catalog Number:

83101-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC052977

GeneID (NCBI):

7133

UNIPROT ID:

P20333

Full Name:

tumor necrosis factor receptor  
superfamily, member 1B

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

230328F5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00777-3, TNFRSF1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83101-4-PBS. Detection antibody: 98010-1-PBS. Standard: Eg0086. Range: 0.625-20 ng/mL.