

For Research Use Only

# Mouse NGAL Recombinant antibody, PBS Only (Capture)



Catalog Number: 83102-6-PBS

## Basic Information

Catalog Number: 83102-6-PBS	GenBank Accession Number: NM_008491.1	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 16819	CloneNo.: 230372D2
Source: Rabbit	UNIPROT ID: P11672	
Isotype: IgG	Full Name: lipocalin 2	
	Calculated MW: 23KD	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Mouse

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00302-3, Mouse NGAL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83102-6-PBS. Detection antibody: 83102-4-PBS. Standard: Eg0572. Range: 0.625-20 ng/mL.