

For Research Use Only

CEACAM8 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83105-1-PBS



Basic Information

Catalog Number: 83105-1-PBS	GenBank Accession Number: BC026263	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 1088	CloneNo.: 230353B7
Source: Rabbit	UNIPROT ID: P31997	
Isotype: IgG	Full Name: carcinoembryonic antigen-related cell adhesion molecule 8	
	Calculated MW: 38 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00293-1, CEACAM8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83105-3-PBS. Detection antibody: 83105-1-PBS. Standard: Eg32133. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00293-2, CEACAM8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83105-4-PBS. Detection antibody: 83105-1-PBS. Standard: Eg32133. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00293-3, CEACAM8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83105-2-PBS. Detection antibody: 83105-1-PBS. Standard: Eg32133. Range: 0.625-20 ng/mL