

For Research Use Only

DDX19A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83106-4-PBS

Basic Information

Catalog Number: 83106-4-PBS	GenBank Accession Number: BC005162	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 55308	CloneNo.: 230232D7
Source: Rabbit	UNIPROT ID: Q9NUU7	
Isotype: IgG	Full Name: DEAD (Asp-Glu-Ala-As) box polypeptide 19A	
Immunogen Catalog Number: AG7362	Calculated MW: 54 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00140-2, DDX19A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83106-3-PBS. Detection antibody: 83106-4-PBS. Standard: Ag7362. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00140-3, DDX19A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83106-5-PBS. Detection antibody: 83106-4-PBS. Standard: Ag7362. Range: 0.625-20 ng/mL.