

For Research Use Only

# DDX19A Recombinant antibody, PBS Only (Capture)



Catalog Number: 83106-5-PBS

## Basic Information

Catalog Number:

83106-5-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7362

GenBank Accession Number:

BC005162

GeneID (NCBI):

55308

UNIPROT ID:

Q9NUU7

Full Name:

DEAD (Asp-Glu-Ala-As) box polypeptide 19A

Calculated MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

230232B8

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00140-3, DDX19A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83106-5-PBS. Detection antibody: 83106-4-PBS. Standard: Ag7362. Range: 0.625-20 ng/mL.