

For Research Use Only

TJP2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83111-2-PBS

Basic Information

Catalog Number: 83111-2-PBS	GenBank Accession Number: BC027592	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9414	CloneNo.: 230474F1
Source: Rabbit	UNIPROT ID: Q9UDY2	
Isotype: IgG	Full Name: tight junction protein 2 (zona occludens 2)	
Immunogen Catalog Number: AG13534	Calculated MW: 134 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00241-1, TJP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83111-2-PBS. Detection antibody: 83111-3-PBS. Standard: Ag13534. Range: 0.625-20 ng/mL