

For Research Use Only

NEDD4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83112-4-PBS

Basic Information

Catalog Number:

83112-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16291

GenBank Accession Number:

BC152452

GeneID (NCBI):

4734

UNIPROT ID:

P46934

Full Name:

neural precursor cell expressed,
developmentally down-regulated 4

Calculated MW:

1319 aa, 149 kDa

Purification Method:

Protein A purification

CloneNo.:

230478D8

Applications

Tested Applications:

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00234-1, Human NEDD4 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83112-2-PBS. Detection antibody: 83112-4-PBS. Standard: Ag16291. Range: 1.25-20 ng/mL.