

For Research Use Only

MMP1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83114-5-PBS

Basic Information

Catalog Number: 83114-5-PBS	GenBank Accession Number: BC013875	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4312	CloneNo.: 230343H6
Source: Rabbit	UNIPROT ID: P03956	
Isotype: IgG	Full Name: matrix metalloproteinase 1 (interstitial collagenase)	
	Calculated MW: 54 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00290-2, MMP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83114-4-PBS. Detection antibody: 83114-5-PBS. Standard: Eg0527. Range: 0.625-20 ng/mL