

For Research Use Only

# CELSR2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83115-2-PBS

## Basic Information

Catalog Number:

83115-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001408

GeneID (NCBI):

1952

UNIPROT ID:

Q9HCU4

Full Name:

cadherin, EGF LAG seven-pass G-type receptor 2 (flamingo homolog, Drosophila)

Calculated MW:

318 kDa

Purification Method:

Protein A purification

CloneNo.:

230319C10

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data