

For Research Use Only

CELSR2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83115-3-PBS

Basic Information

Catalog Number:

83115-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001408

GeneID (NCBI):

1952

UNIPROT ID:

Q9HCU4

Full Name:

cadherin, EGF LAG seven-pass G-type receptor 2 (flamingo homolog, Drosophila)

Calculated MW:

318 kDa

Purification Method:

Protein A purification

CloneNo.:

230319D10

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00165-1, Human CELSR2
Recombinant Matched Antibody Pair - BSA and
Azide Free measured by Cytometric bead array.
Capture antibody: 83115-4-PBS. Detection
antibody: 83115-3-PBS. Standard:null. Range: null.