

For Research Use Only

# INO80 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83117-5-PBS



## Basic Information

Catalog Number:

83117-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC146785

GeneID (NCBI):

54617

UNIPROT ID:

Q9ULG1

Full Name:

INO80 homolog (S. cerevisiae)

Calculated MW:

1556 aa, 177 kDa

Purification Method:

Protein A purification

CloneNo.:

230440D4

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00159-3, Human INO80 Recombinant Matched Antibody Pair - PBS only. 83117-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY00250. 83117-5-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean INO80 concentration was determined to be 4.4 ng/mL in Jurkat cell extract based on a 2.3 mg/mL extract load.