

For Research Use Only

GCNT1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83121-2-PBS

Basic Information

Catalog Number: 83121-2-PBS	GenBank Accession Number: BC074885	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2650	CloneNo.: 230447G12
Source: Rabbit	UNIPROT ID: Q02742	
Isotype: IgG	Full Name: glucosaminyl (N-acetyl) transferase 1, core 2 (beta-1,6-N- acetylglucosaminyltransferase)	
Immunogen Catalog Number: AG12264	Calculated MW: 428 aa, 50 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00296-1, GCNT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83121-2-PBS. Detection antibody: 83121-3-PBS. Standard: Ag12264. Range: 0.625-20 ng/mL